



CERKAMED

MEDICAL COMPANY

MADE IN POLAND

PRODUCT CATALOG



CERKAMED.COM

...because details are all-important



CERKAMED

MEDICAL COMPANY

We are
Polish producer
of a dental materials



25

Years of experience

KEY TO SUCCESS



High quality products



Wide range of products



Products equipped with accessories

Fast delivery of goods



Rapid order processing

Customized special offers



Excellent shipment protection

Low level of complaints



Stalowa Wola

WE GUARANTEE HIGH QUALITY

and safety of our products through a certified
Quality Management System for Medical Devices, which
meets the requirements of the **EN ISO 13485** standard.



Product sales to over
+100
 countries

all over the world, i.a. in Europe,
 Asia, Africa and South America

We actively
 participate in the
 largest dental trade
 fairs in the **WORLD.**



The reliability of our services is confirmed by numerous awards and certificates. Each subsequent distinction fills us with pride and it remains a guarantee of compliance with restrictive registration procedures around the world.



NOVELTIES 2026

LUMINESCENCE WHITE

Professional teeth whitening gel
(Melon, Mint, Neutral, Strawberry)

16% carbamide peroxide

Key features of LUMINESCENCE WHITE

- ✓ 16% carbamide peroxide – releases approx. 5.8% hydrogen peroxide
- ✓ Contains potassium nitrate and sodium fluoride (PF)
- ✓ For use with custom-made or universal thermoplastic trays
- ✓ Recommended application: 4–6 hours once daily
- ✓ Recommended treatment period: approx. 10 days
- ✓ 1.2 ml syringe provides approx. 2–3 applications
- ✓ Available in 4 variants: Neutral, Melon, Mint and Strawberry

2x1,2 ml



2x1,2 ml
+ tray case



cerkamed.com



POLISH BOOST

Professional polishing paste
(Sensitive Strawberry, Regular Mint,
Intensive Neutral)



ACT combines:

- ✓ Precisely selected ingredients
- ✓ Structured mode of action
- ✓ Controlled and predictable polishing
- ✓ Precisely defined abrasives



5

Key features of POLISH BOOST


- ✓ Professional cleaning and polishing paste for use in dental practices
- ✓ Medium abrasiveness – RDA 40
- ✓ Developed with ACT – Active Crystal Technology
- ✓ Available in a 50 g tube



RUBBER-DAM LIQUID

Pearl White

Light-curing Rubber-Dam Liquid



1,2 ml


- ✓ Gingival protection during the treatment
- ✓ Sealing rubber dams
- ✓ New colour – Pearl White



TOP FLEX

Dental irrigation tips

20 pcs



0.3 mm (30G)



- ✓ Flexible cannula
- ✓ 27 mm working length
- ✓ Radial irrigation
- ✓ Depth markings

CERKAMED TOOLS

New tools, new possibilities.
28 variants, available soon,
stay tuned!

Restorative Dentistry

NEW IN:
28
TOOLS

Endodontics



Comfort

Corrosion resistance

Durability

Perfect sharpness

AVAILABLE SOON



FRESH DENT SMART

Ready to use mouthwash (Mint, Original, Kids, Melon Mint)

7

Key features of FRESH DENT SMART

- ✓ Developed with ACT – Active Crystal Technology
- ✓ Ready to use, handy 100 ml bottle – perfect for your bag, work or travel
- ✓ Small size, everyday convenience – easy to keep close at hand
- ✓ Quick 30-second rinse for a fresh and clean feeling
- ✓ Available in 4 variants: Mint, Original with Chlorhexidine and Melon Mint, Melon - dedicated for kids



100 ml



SOFT GOGGLES

Soft protective eyewear for a patient

10 pcs



SPIT SPOT

Disposable dental cuspidor liners

50 pcs



PROTECTIVE PADS GHOST

For a patient

30 pcs



MOUTH OPENER

(S, M, L)

2 pcs set



Silicone

50%
MORE PRODUCTS




CHLORAXID

sodium hypochlorite (2% or 5,25% of active chlorine), purified water
compliant with Ph. Eur.
Liquid for root canals rinsing

Available packages:

- CHLORAXID 2% 200 g
- CHLORAXID 2% 400 g
- CHLORAXID 5,25% 200 g
- CHLORAXID 5,25% 400 g



ENDO-SOLUTION

EDTA disodium salt (17%)
LIQUID for root canals widening

Available packages:

- Bottle 200 g



CITRIC ACID 40%

40% citric acid solution
LIQUID for root canals rinsing

Available packages:

- Bottle 200 g
- Bottle 400 g






ENDO-PREP CREAM

EDTA disodium salt (15%), urea peroxide 10%
CREAM for root canals preparation

Available packages:

- Syringe 2 ml
- Syringe 5 ml
- Syringe 10 ml
- MEGA PACK 3 syringes x10 ml



ENDO-PREP CREAM



ENDO-PREP GEL

EDTA disodium salt 17%
GEL for root canals preparation

Available packages:

- Syringe 5 ml
- Syringe 10 ml





ENDO-TOP

Endo irrigation needles

Available variants:

- Box 100 pcs



CE
2409

CAVIPACK

Zinc Oxide, Zinc Sulfate, Calcium Sulfate, excipients
Temporary filling

Available packages:

■ Package 3x10 g

- ready to use, zinc oxide and zinc sulfate paste for temporary filling
- it does not contain eugenol
- binds chemically when exposed to moisture
- easy to remove
- radiopaque
- pleasant, rose scent
- economy package **3x10g**



CE
2409

CANAL DETECTOR

Methylene blue
Root canals orifices indicator

Available packages:

■ Syringe 2 ml



CE
2274

CALCIPAST

Calcium hydroxide
Calcium hydroxide paste for temporary root canal filling

Available packages:

- Syringe 2,1 g
- MEGA PACK 4 syringes x2,1 g

before

after



CE
2409

RED DETECTOR

Rhodamine B
Caries detector

Available packages:

■ Syringe 2 ml



CE
2409



ZINC OXIDE FAST

Zinc oxide, excipients
Fast-setting root canal repair material in combination with eugenol

Available packages:

□ Container 50 g



ZINC OXIDE

Zinc oxide
Material for making pastes in combination with eugenol

Available packages:

□ Container 50 g

CE
2409



EUGENOL

Eugenol oil 100%
Oil for making pastes

Available packages:

- Bottle 10 ml
- Bottle 20 ml




CE
2409

MTA+

Material for filling
and rebuilding of root canals

Available packages:

- MINI 3 x 0,14 g MTA+ powder,
1 ml liquid
- MAXI 10 x 0,14 g MTA+ powder,
1 ml liquid



MTA



CE

MTA+ APPLICATOR

Reusable dental tool for handling
and application of prepared
measure of MTA+ material

Possible variants:

- 0,8 mm
- 1,2 mm



CE

TOOLS FOR HANDLING AND APPLICATION OF MTA+

Reusable dental tools for handling
and application of prepared
measure of MTA+ material



MTA SMART SYSTEM TRAINING SET

Endodontic training set
for application MTA

SET INCLUDES:

SILICONE TRAINING MASS

TRAINING TEETH 5 PCS

SILICONE TRAINING MASS STAND
ALLOWS 3D ROTATION

TRAINING
MTA POWDER 5 g

TRAINING
MTA LIQUID 3 ml

MEASURING CUP
FOR MTA POWDER

SLIDE

PIPETTE FOR
MTA LIQUID APPLICATION



CE
2409

BIO MTA+

Biomaterial for filling
and rebuilding of root canals

Available packages:

- MINI 3 x 0,14 g BIO MTA+
powder, 1 ml liquid
- MAXI 10 x 0,14 g BIO MTA+
powder, 2 x 1 ml liquid

TRAINING TEETH



2633-101-PCL

Maxillary first molar
3 canals, 3 roots
Perforation of the
bottom of the chamber
floor, lateral wall
perforation of the palatal
root



2111-1A2-PCL

Maxillary central incisor
1 canal, 1 root
Incomplete root
development



2633-102-PCL

Maxillary first molar,
3 canals, 3 roots,
perforation of the pulp
chamber floor, external
apex resorption



3623-101-PCL

Mandibular first molar
3 canals, 2 roots
perforation of the
pulp chamber floor



3111-1A1-PCL

Mandibular central
incisor 1 canal, 1
root, anterior root
wall perforation

Also sell
separately!

DRSK

Modified DRSK RCT tooth model

5 variants available

CE
2409

SYNTEX

Root canals filling
and sealing material

Available packages:

- Syringe 10 g
+mixing tip
+mixing papers




CE
2409

GUTTA PERCHA POINTS

Hand rolled gutta percha points

Available packages:

- Package 2% 120 pcs marked
- Package 4%, 6% 60 pcs marked
- Package ProTaper 60 pcs marked
- 2% ROUND PACKAGE 240 pcs marked
- 4%, 6% ROUND PACKAGE 240 pcs marked



CE
2409

GUTTA PERCHA PLUGGER

Reusable dental tool for vertical
condensation of gutta-percha

Available packages:

- Single package
- Set of 4 sizes



CE
2409

ABSORBENT PAPER POINTS

Hand rolled absorbent paper points

Available packages:

- Package 2% 200 pcs marked
- Package 4%, 6% 100 pcs marked
- Package ProTaper 100 pcs marked
- 2% ROUND PACKAGE 400 pcs marked
- 4%, 6% ROUND PACKAGE 240 pcs marked



CE
2409

EUCALYPTOL

Eucalyptus oil 100%
Oil for dissolving gutta-percha

Available packages:

- Bottle 10 ml
pipettes for application 3 pcs



CE
2409

ORANGE GUTTANE

Orange oil 100%
Oil for dissolving gutta-percha

Available packages:

- Bottle 10 ml
pipettes for application 3 pcs



ENDO aspirator^{PRO}

for removal liquids from the root canals

Available variants:

- ENDO-ASPIRATOR PRO
- ENDO-ASPIRATOR PRO Golden Version


PACKAGE INCLUDES:

- HANDLE - 1 pc
- CONNECTOR - 3 pcs
- LUER LOCK HOSE - 5 pcs
- MAGNET - 5 pcs
- RAMROD FOR CLEANING THE ENDO-ASPIRATOR PRO - 1 pc



ENDO-ASPIRATOR

for removal liquids from the root canal
by aspirating



ADDITIONAL ACCESSORIES

PACKAGE INCLUDES:

- LUER LOCK HOSE - 5 pcs
- CONNECTOR - 1 pc



btr^{Pen} ^{CE}

Dental tool for removing broken canal instruments

Set includes:

- BTR PEN – 1 pc stainless steel tool
- BTR CALIBRATOR – 1 pc
- BTR TIP Ø 0,3 – 5 pcs
- BTR SPRING – 2 pcs soft and hard (third middle spring fixed on the device)
- Silicone BTR O-RING – 3 pcs



precise grip of the broken instrument

life-long warranty

single-tool system for every type of root canal

btr-pen.com



THREE TIPS ONE GOAL

Ø 0,3 mm

wire diameter 0.076 mm



Ø 0,4 mm

wire diameter 0.085 mm



Ø 0,5 mm


wire diameter 0.1 mm



Changeable tips available separately, packed in 5 pcs



btr^{Pen}
OFFICIAL GROUP




L I M I T E D

Golden edition

^{CE}

Set includes:

- BTR PEN Limited Golden Edition – 1 pc stainless steel tool
- BTR CALIBRATOR – 1 pc
- GOLDEN BTR TIP – 0,3 mm – 20 pcs
- BTR TIP – 0,4 mm – 10 pcs
- BTR TIP – 0,5 mm – 10 pcs
- BTR SPRING – 2 pcs soft and hard (third middle spring fixed on the device)
- Silicone BTR O-RING – 3 pcs



btr-pen.com

Join the world's endodontic elite
with limited edition of BTR PEN

BTR DENTAL TRAINING KIT

Ideal kit for training the removal of broken endodontic files from root canals

1. DIGITAL MICROSCOPE ADSM201

- MINI HDMI input
- MINI USB port
- MICRO SD card input
- light regulation
- LED lightning

2. MICROSCOPE'S INSTRUCTION FOR USE

3. TWO TRAINING TEETH

with pieces of broken endodontic files

HARD 3 deeply placed, blocked pieces of files (requiring processing before removal)

EASY 3 loosely placed pieces of files

4. SILICONE TRAINING MASS

- transparent and backlit
- 8 pieces of broken endodontic files
- a space to place the training tooth

5. SILICONE MASS STAND

- allows 3D rotation

6. BTR Pen WORK METHODOLOGY

- work with BTR Pen step by step

7. ANTI-SLIP MAT

- stability during the work with the kit



BTR DENTAL TRAINING GUM

1. SILICONE TRAINING MASS

with pieces of broken endodontic files

2. SILICONE TRAINING MASS STAND

- allows 3D rotation

3. TWO TRAINING TEETH

with pieces of broken endodontic files (easy/hard)

4. ANTI-SLIP MAT

DRSK
Modified DRSK RCT tooth model



Join our group



btr-pen.com



facebook.com/groups/btrpen



BLUE ETCH FLOW

O-phosphoric acid 36%

Flowable dental etchant

Available packages:

- Syringe 2 ml
- Syringe 10 ml
- Syringe 50 ml
- MONSTER PACK 10 syringes x 2 ml





BLUE ETCH

O-phosphoric acid 36%

Dental etchant

Available packages:

- Syringe 2 ml
- Syringe 10 ml
- Syringe 50 ml
- MEGA PACK 3 syringes x 10 ml
- MONSTER PACK 10 syringes x 2 ml



YELLOW PORCELAIN ETCH

Hydrofluoric acid 9,5%

Porcelain etch

Available packages:

- Syringe 2 ml



SILAN

Methacrylate silane

Available packages:

- Syringe 2 ml





RAINBOW FLOW

Methacrylate resins
Light-curing flow composite

Available packages:

■ syringe 1 g (chosen colour)
+ applicators

- orange
- blue
- pink
- yellow
- green
- purple



3D
Stickers



RAINBOW FLOW dentina

Methacrylate resins
Light-curing flow composite

Available packages:

■ Syringe 1 g + applicators

Application:

▼ Pit and fissure sealing



▼ Filling of permanent teeth and
milk teeth or protection of root
canal orifices



A2

filler content 62%



FLUOR DEFENDER

0,1% fluorosilan (1600 ppm F)
HEMA

Protective enamel with
fluoride for children

Available packages:

- MINI 3 vials x 1 ml
- + set of applicators
- MAXI 10 vials x 1 ml
- + set of applicators



MARKETING MATERIALS FOR PATIENT

Available packages:

1. 3D STICKERS 10 pcs, 50 pcs
 2. PUPPET 3 pcs
 3. BOXES FOR MILK TEETH 5 pcs
 4. DIPLOMA „BRAVE PATIENT“ mix 20 pcs
 5. GUESSING GAMES FOR A BRAVE PATIENT 20 pcs
 6. TEETH BRUSHING CALENDAR for kids
- Set included: teeth brushing calendar, diploma, stickers, magnet



ALUSTAT FOAM

Aluminium chloride 20%
FOAM for staunching the bleeding

Available packages:

- Syringe 0,8 g
- MEGA PACK 4 syringes x 0,8 g



CE 2409

MEGA PACK

ALUSTAT FOAM
Foam for staunching bleeding
TO BE USED BY DENTIST ONLY

aluminium chloride 20%
CE 2409



ALUSTAT

Aluminium chloride 25%
LIQUID for staunching the bleeding

Available packages:

- Bottle 10 g



ALUSTAT GEL

Aluminium chloride 25%
GEL for staunching the bleeding

Available packages:

- Syringe 5 ml
- Syringe 10 ml

TCA DENT

Trichloroacetic acid 35%
Preparation for local cauterisation of the mucous membranes

Available packages:

- Bottle 10 ml





BEST-CORD NANO

Knitted retraction cords,
254 cm lenght

Polyester/Polyamide yarn conjugated.
Aluminium chloride impregnated:
000.0/000/00/0: $0,5 \pm 0,3$ mg/inch
1/2: $1,0 \pm 0,4$ mg/inch

6 SIZES:

000.0
000
00
0
1
2

CE 2409



3 TYPES



RETRACTION CORD PACKERS

Tools for packing retraction cords in the gingival pockets

Possible variants:
RETRACTION CORD PACKERS

- UNIVERSAL 1,5 mm
- ROUND 1,5 mm
- ROUND 3 mm



GLUCOSITE

aqua, chlorhexidine digluconate, alcohol,
hydrogen peroxide

LIQUID for gingival pockets

Available packages:

- Syringe 2 ml
- Bottle 50 ml



GLUCOSITE GEL

aqua, glycerine, chlorhexidine digluconate,
methylcellulose, alcohol, hydrogen peroxide

GEL for gingival pockets

Available packages:

- Syringe 2 ml
- Bottle 50 ml



PULP SPRAY

SPRAY for pulp vitality testing

Available packages:
Spray 200 ml

Aroma:

- MINT
- ORANGE



INOX

propanotriol (glycerol), Ph Eur
purified water, gel base

**Gel promoting complete
polymerisation of dental
composites**

Available packages:

- syringe 2 ml
- MEGA PACK
4 syringes x 2 ml

Each item of
Rubber dam Set
available

separately!

RUBBER-DAM SET

Tools for rubber dam Start Pack



Includes:

- RUBBER-DAM PUNCH
- CLAMPS FORCEPS
- RUBBER-DAM FRAMES (3 sizes)
- THE SET OF CLAMPS
- BASE FOR CLAMPS

**12 types
of clamps**

Also available in satin version

3

13A

208

5

W8A

12A

2A

210

212

27

29

00

Flat jaws for small lower molars (No 3).
Festooned wings for large upper molars.
Wings fit the gum curvature (No 5).
Festooned jaws for molar's roots (No W8A).
Serrated jaws for lower right and upper left molars (No 12A).
Serrated jaws for lower left and upper right molars (No 13A).
Small and large jaw for anteriors (No 210).
Flat jaws for anteriors (No 212).
Rounded jaws for premolars (No 27).
Rounded and flat jaws for premolars (No 29).
Small jaws and high bow for anteriors and premolars (No 00).
Wings, flat and large jaws for premolars (No 2A).
Festooned wings, flat jaws for large premolars (No 208).

RUBBER-DAM



Available package: 36 pcs

Available variants:



medium 0,20 mm

**mint
BLACK**

medium 0,18 mm
or heavy 0,22 mm

**plain
BLUE**

medium 0,18 mm
or heavy 0,22 mm

**mint
GREEN**

NON LATEX DENTAL DAM



Available package: 36 pcs

Available variants:



NEW

0,18 - 0,23 mm

MEDIUM

0,23 - 0,28 mm

HEAVY



RUBBER-DAM LIQUID

urethane resins, inorganic fillers

Light-Cured RUBBER-DAM liquid

Available packages:

- Syringe 1,2 ml
- MEGA PACK 4 syringes x 1,2 ml



RUBBER-DAM TEMPLATE

Enables easy hole making for teeth in dental rubber dam

Size:

■ 6x6 – 150x150 mm

RUBBER-DAM NAPKINS

Available packages:

- Set of 36 pcs



RUBBER-DAM CORD

Dental dam stabilizing cord

Available packages:

- RUBBER-DAM CORD fine
- RUBBER-DAM CORD medium
- RUBBER-DAM CORD thick

Package containing 1,2 m of rubber dam cord

FINE



MEDIUM



THICK





ISOPROPYL ALCOHOL

Available packages:

- Bottle 200 g



IODOFORM

Available packages:

- Container 30 g



METHYLENE CHLORIDE

Available packages:

- Bottle 10 ml


DENT FRESH

Available packages:

- Start Pack
concentrate
in syringe 50 ml
+ bottle
+ dispenser
- concentrate
in syringe 50 ml

Two variants available:

- MINT
- ORIGINAL



CALCIUM HYDROXIDE

Calcium hydroxide pure p.a.

Available packages:

- Container 10 g

Perfect picture!

15
pcs

Anti-reflex

DENTAL CONTRASTOR

universal
single use
flexible



SILICONE RINGS

Silicone rings for marking
of instruments in stomatology

Available packages:

- Box 20 pcs
- Box 100 pcs



Available colours:

- black
- purple
- blue
- red
- green
- yellow
- white

PROTECTIVE VISORS



Available colours:

- green
- blue
- pink
- beige
- white
- transparent

COMFORT LIGHT
with set of foils



STANDARD / RELAX
with set of foils



+Set of 10 foils!



PROTECTIVE VISORS ATTACHED TO GLASSES





Available package:

- The set contains 2 pieces of reusable foils that can be washed



Comfort of work!




ANTISTATIC FOILS



SET of antistatic interchangeable foils

Available packages:

- 4 x COMFORT RELAX (0,175 mm)
- 5 x COMFORT STANDARD (0,125 mm)
- 4 x COMFORT LIGHT (0,175 mm)



EASY 2

PROTECTIVE SHIELDS FOR A PATIENT



Packing includes:

- Shield frame with plastic hooks + set of exchangeable foils - 8 pcs

Available colours:

- green
- blue
- pink
- beige
- white
- transparent



EASY 3

PROTECTIVE PADS FOR A PATIENT



Packing includes:

- 20 pcs OF DISPOSABLE PADS

one size



GLASSES FOR POLYMERIZATION



MINI

ENDO-STATION GO



Line of endodontic organizers


Available colors*:

- GREEN
- BLUE
- TRANSPARENT


* ENDO-STATION GO GUTTA BOX available only TRANSPARENT



GUTTA BOX 12



GUTTA BOX 2x6



GUTTA BOX 2x3

Przedsiębiorstwo Produkcyjno-Handlowe

CERKAMED

WOJCIECH PAWŁOWSKI



ul. Kwiatkowskiego 1
37-450 Stalowa Wola
Poland



+48 15 840 08 09
+48 15 842 35 85, int. 809
+48 785 665 505



export@cerkamed.pl



www.cerkamed.com

EN version



Find
us on
Facebook



Follow
us on
Instagram

btrPen
OFFICIAL GROUP